



Polycarbonate Erlenmeyer Culture Flasks

Ideal replacement for glass and other Erlenmeyer/shaker flasks.

- Patented [DuoCAP®](#) allows sterile air exchange or a leak-proof seal.
- Minimize oxidation or pH shifts.
- Perform aerobic/anaerobic culture.
- Available with [septum cap](#) for pathogen containment.
- Re-autoclavable, reusable.
- Individually Packed in EZ-Tear bags. Sterile to SAL 10⁻⁶
- Available with a flat or a baffled base in desired packaging options.
- Choose baffled base flasks when increased levels of oxygenation add mixing are required



Leakproof



Optically Clear



Shatterproof



Autoclavable



Gamma Radiated



BioPioneer Inc.

8540 production Ave. Suite A San Diego CA 92121

Tel: 858-228-6949, Fax: 858-225-0288

<http://www.BioPioneerinc.com> Info@BioPioneerinc.com