

Product Description

100bp DNA ladder

MDL-100, 100 applications

100 bp - 1.5kb DNA Ladder , ready to use

11 DNA fragments: 100, 200, 300, 400, 500x3, 600, 700, 800, 900, 1000x2 and 1500bp


For 100 applications

Concentration: 800ng / 8 μ l

Recommended loading volume: 5 μ l - 8 μ l

Storage Buffer: 10mM Tris-HCl, pH8.0, 1mM EDTA, and 1 X Loading Buffer (5% Glycerol, 0.02% bromophenol blue, 0.02% xylene cyanol)

Storage -20°C



BioPioneer Inc.

8540 Production Ave., Suite A, San Diego, CA 92121

Tel: 858-689-6988 FAX: 858-225-0288

Email: info@Biopioneerinc.com Web: <http://www.Biopioneerinc.com>